



Classical beauty is a modern conquest.
Three times more beautiful with
Oxygen, Radiofrequency and Ultrasounds



TRIO2

Le bellezza é un fatto scientifico:

- Ultrasuoni
- Radiofrequenza
- Trattamenti all'ossigeno per viso e corpo.

Beauty is a matter of science:

- Ultrasounds
- Radio-frequency
- Oxygen facials and body treatments.





Firming inside, Reshaping and Revitalizing outside. What else would you like from a face & body system?



- TriO2 stimulates and firms deeply the tissue by means of the reparative processes started by **radiofrequency**
- It reshapes through the powerful anti-cellulite action of **ultrasounds**
- It oxygenates and revitalizes the skin, driving deep into the derma 93% pure **oxygen** together with natural products based on wheat grass and other studied fresh plants' extracts

TriO2 comes from the twenty-year-old Exar experience in the “oxygen for beauty and wellness” field. TriO2 opens the road for interesting contaminations, putting together three precious functions: Radiofrequency, Ultrasounds and Oxygen.

1. **Monopolar Radiofrequency at 460 KHz – for face & body firming.** The R.F. function let us stimulate deeply the tissue – thanks to its soft thermal action it produces the denaturation of collagen producing a lifting effect and it increases the production of new collagen fibers stimulating the natural skin processes. For visible long lasting skin tighteninging results.
2. **Ultrasounds at 1 MHz –** TriO2 sends sound waves (at 3 million vibrations per second) very deeply into the skin in order to reduce cellulite, enhance local blood flow, nourish all the skin layers, improve tissue tone and carry into the skin active principles (sonophoresis). Note - considered that the penetration becomes lower as the frequency increases, the high-frequency Ultrasounds of TriO2 are perfect to cope successfully with cellulite.
3. **Oxygen facials and Oxygen body treatments** thanks to the Oxygen Spray function that let us gently spray on the skin pure oxygen together with natural light products – all highly nourishing and rich in antioxidants. For effective not-invasive beauty treatments – able to ensure fast concrete results on face and body.



Technical SPEC

Dimensions	52cm x 52cm x 95/107cm (H min/max)	
Weight	40 kg	
Voltage	220 V \pm 15%	
Power Consumption	min. 100 Watt (U.S. function), 200 Watt (R.F. function), 350 Watt (Oxygen function)	
External protection	2 fuses 6,3 A-250 V	
Classification	Class 1 – Type BF (Ref: CEI 62-39)	
Temperature	10°~40°	
OXYGEN F unction Mode oxygen production: PSA system (Pressure Swing Adsorption)		
Oxygen spray	Oxygen Purity	90% \pm 3%
	Oxygen Flow Rate / Minute	6L/Min \pm 0,5
	Handle	OxyGun Sprayer Nozzle
RADIOFREQUENCY F unction Emission mode: monopolar capacitive		
Monopolar capacitive Radiofrequency	Power input	230 V +/- 10% 50-60 Hz
	Max power absorbed	200 Watt
	R.F. Frequency	470 KHz +/- 10%
	Max power at the exit	170 Watt +/- 10%
	Independent Exit	1
	R.F. handle	1
	Electrodes (black lacquered in order to be electrically isolated) that can be connect to the R.F. handle	6 (from diameter 2,5 cm to 7,5 cm)
	Mass (grounding)	1 steel plate
	Foot pedal to adjust R.F.	1
	Programs	Face and Body, customizable
ULTRA SOUNDS Function High Frequency Ultrasounds		
High Frequency Ultrasounds	Power input	220V – 50-60Hz
	Max power absorbed	100 Watts
	U.S. Frequency	1 MHz
	Max power at the exit	30 Vpp max, 2 mA max
	U.S. Handle	1
	Programs	8 body (4 against cellulite, 4 against adipose tissue)